

**NEW**



**AARNA SYSTEMS**

# Body Composition Analyser



Jawon Medical Co., Ltd has earned CE (Europe), FDA (USA) and KFDA (Korea). Its Quality and achievement have been approved by these certificate.





The technology of our company makes the professional body composition analyzer **light weight** and **performs the whole function.**

**SENSATIONAL DESIGN**

**RELIABILITY**

**CONVENIENCE**



## VISUAL COMMUNICATING FOR HEALTH

Body Composition Analysis  
Assessment of Weight Control  
Abdominal Analysis

Body Type  
Segmental Assessment  
Control Guide  
Body Composition Change

 **JAWON MEDICAL**  
Our company name will be changed to JANEX MEDICAL.



# ULTRA LIGHT BODY COMPOSITION ANALYZER FOR HOSPITALS

## SUPERIOR FUNCTION WITH COMPACT SIZE BODY COMPOSITION ANALYZER

Main Body  
Approx.  
**10kg**



Designed by special structure



### CARRY BAG Type A

ioi can be transferred easily and used everywhere owing to ioi exclusive bag.

Size of carry bag  
(approx. 815 × 215 × 410)mm  
Weight of carry bag : approx.1kg  
(option)

Weight of carry bag + main body

**Approx. 11kg**



### CARRY BAG Type B


The exclusive plastic bag has wheels and string in order to move easily and safety.

Size of carry bag  
(approx. 866 × 295 × 567)mm  
Weight of carry bag : approx.6kg  
(option)

Weight of carry bag + main body

**Approx. 16kg**





7 inch wide color LCD

The technology of our company makes the professional body fat measuring device **light weight** and **performs the whole function**.

01. Carry bag for safe mobility.(option)
02. Designed by human engineering
03. It has been designed so that it will not be damaged when dropped at a height less than 50cm. (carry bag)
04. Ultra light, approx 10kg, body composition analyzer
05. You can use both A4(basic) and thermal(option) printers.
06. 7inch wide color screen
07. Non-bending weighing scale frame
08. Easy installation, 3-5 minutes, can be done by anyone.
09. Mobile friendly, it can be located anywhere.
10. Easy measurement process through Voice Guide.
11. Saving of measured data by USB Memory Stick.(option)
12. It is convenient when going on a business trip or being serviced as it is mobile friendly.

• Basic offer



Management program for body composition data



PRINTER & CART

Model or specification of printer & cart can be subject to change according to market demand.

• ioi353 Option



Professional counseling program for practical use

Easy Body Plus



USB M

### 01. 5 factors

Jawon Medical Co., Ltd has been using five factors in the calculation of the body composition since the development stage in 1996. It has been rewarding to see the technology progress to its current level.

#### What is Bioelectrical Impedance Analysis? (BIA)


In BIA, the difference in electrical conductivity according to biological characteristics of the tissue has some limitations as follows.

#### The limitation of BIA:

The human body is represented as a simple cylindrical shape determined by height and weight. It assumes that body composition is homogenous and evenly distributed. It does not consider individual differences and variation of Body Composition. It ignores any changes in the environment (temperature), body heat, and stress.

However, in reality, the human body is different with the assumptions mentioned above.

In order to use BIA to analyze the "actual" human body accurately, 5 factors must be considered in the calculation.



Jawon medical uses 5 factors correctly to analyze the whole body and abdominal fat.

Using 5 factors, the accuracy of final analysis is greatly improved.

### 02. Tetra Polar Method


We, Jawon Medical Co., LTD use the Tetra Polar Method to achieve the greatest accuracy in analyzing body composition.

Tetra Polar Method uses separate electrode;


Current Electrode which sends out an electric current to the human body and Voltage Electrode which detects impedance.

It minimizes the contact resistance and increases the accuracy.

In addition, it takes a measurement at the fixed location of foot electrode (ankle electrode) and hand so that it reduces measuring errors and increases reproducibility of results.



### 03. Multi Frequency



**FOR YOUR HEALTHY & BEAUTIFUL LIFE**  
ioi 353 can analyze correctly human body's mystery that cannot see externally.

**ioi 353** BODY COMPOSITION ANALYSIS

**JAWON MEDICAL**  
www.jawon.com


Name / ID Lauren 000000001  
Date 2008-06-08 16:25:11  
Height 155.4 cm Age 49 yrs  
Weight 69.0 kg Gender Female

**Body Composition** [kg]

Weight / Over 69.0 [47.8~58.4]	Std.wt. 53.1	Body Fat 25.2	
L.B.M. / Over 43.8 [37.2~42.5]	Mineral 3.5		Body Fat 25.2
S.L.M. / Over 40.3 [34.2~39.3]	T.B.W. / Over 31.5 [26.7~30.6]		Body Fat / Over 25.2 [10.6~ 15.9]
Protein / Over 8.8 [7.4~8.5]		Mineral / Over 3.5 [2.9~3.1]	

Std.wt. : Standard weight L.B.M. : Lean Body Mass S.L.M. : Soft Lean Mass T.B.W. : Total Body Water  
Mineral is conservative estimate.  
The assessment of Under, Optimal and Over is decided by standard weight on Body Composition table.

**Assessment of Weight Control**



B.M.I. : Body Mass Index P.B.F. : Percent Body Fat

**Abdominal Analysis**



A.C. : Abdominal Circumference  
A.C. is an estimated value in case of measure the navel circumference.  
Accuracy of abdominal analysis may be decreased in case of idiowascularity.

**Control Guide** [kg]


	Measured data	Control	Goal to control
Weight	69.0	+15.9	Target to control +11.9
M.B.F.	25.2	+11.9	Control/week 0.5
S.L.M.	40.3	+3.6	Duration of control 24 week
B.M.R.	1248 kcal	T.E.E.	1 922 kcal
A.M.B.	53 yrs Impedance	445 Ω	

M.B.F. : Mass of Body Fat B.M.R. : Basal Metabolic Rate T.E.E. : Total Energy Expenditure A.M.B. : Age Matched of Body  
Age Matched of Body is reference value.  
Control guide and calorie prescription are proposed value for your body type.

**Body Type**



**Segmental Assessment**



Item	Lt.Arm	Rt.Arm	Trunk	Lt.Leg	Rt.Leg
M.B.F. (kg)	1.83	1.69	13.76	4.10	3.82
S.L.M. (kg)	2.58	2.60	20.20	7.45	7.47


Segmental Assessment is reference value.

**Blood Pressure**

Systolic mmHg	137	Diastolic mmHg	86
Pulse bpm	86		

P.R.P. : Pressure Rate Product

You need to control 550 kcal from T.E.E. 1922 kcal.



**Body Composition Change**

	Date	Weight	M.B.F.	S.L.M.
Previous	08.05.25	69.6	25.5	40.0
Present	08.06.08	69.0	25.2	40.3

Rev. A.0(08.07.29)


## Environmental parameters

Setup function	Contents
ID usage	It is selected whether ID is used for subjects or not.
Scale offset	Compensating measured value of weight scale
Clothes	Compensating the weight of clothes worn
Print position	Adjusting print position to fit to the pre-formatted result sheet in the direction of up/down and left/right
Date/Time	Setting current date and time
Logo	Printing hospital name, address, contact information, and logo

## Specification

Model	ioi 353
Measuring method	BIA via tetra-polar electrode method using 8 touch electrodes
Frequency range	5, 50, 250 kHz
Measuring site	Whole body and Segmental measurement (right arm, right leg, left arm, left leg, and trunk)
Result item	Protein, Mineral, Total Body water, Mass of Body Fat, Soft Lean Mass, Lean Body Mass, weight, standard weight, B.M.I., Percent Body Fat, Age Matched of Body, Basal Metabolic Rate, Total Energy Expenditure, Body type, 5 body parts (right arm, right leg, left arm, left leg, and trunk) Soft Lean Mass/Mass of Body Fat and assessment, Body Composition Change, Control guide (weight/Mass of Body Fat/Soft Lean Mass Control, Goal to control, Control/week, Duration of control, Diet prescription calorie, Exercise prescription calorie), Visceral Fat Area, Visceral Fat Level, Abdominal Circumference, W.H.R., Impedance, Blood pressure (when connected with blood pressure monitor of our company)
Current	Within 280μA
Power consumption	60VA
Power supply	Input: AC 100~240 V, 50/60 Hz, Output: DC 12 V, 5 A adapter
Display	7 inch wide color LCD
Input device	Key pad, PC remote control
Transmitting device	USB port, RS-232 port (Wireless)
Printing device	USB port, Thermal Printer(option)
Dimension	400 x 735 x 890mm (W x D x H, ± 10mm)
Weight	About 10kg (main unit)
Measuring range	100~950 Ω
Measuring time	Within 1 minute
Input height	100~200cm
Measuring weight	10~200kg
Input age	5~89 years old
Operation ambient	Temperature 10~40°C, Humidity 30~75% (non-condensing)
Storage ambient	Temperature -20~60°C, Humidity Less than 95% (non-condensing)

※ For purposes of improvement, specifications and design are subject to change without notice.



Certification of Jawon



**JAWON MEDICAL**  
www.jawon.com

JAWON MEDICAL CO.,LTD(The Head Office)  
FACTORY / 1208-12, Shinsang-Ri, Jinryang-Eup, Kyungsan-City, Kyungsang Bukdo, 712-837, KOREA  
TEL : 82-53-856-0993 FAX : 82-53-856-0995 E-mail : jawon@jawon.com

FACTORY II / 1186-8, Shinsang-Ri, Jinryang-Eup, Kyungsan-City, Kyungsang Bukdo, 712-837, KOREA  
Our company name will be changed to JANEX MEDICAL.

IOI 353



**AARNA SYSTEMS**

Established in 2008, we Aarna Systems are renowned Importer, Trader, Supplier and Distributor (in India) of Ultrasound Bone Densitometer, Body Composition Analyzer and Portable BMD Machine.

We are importer & distributor (in India) of Ultrasound Bone Densitometer CM-200 of Furuno Electric Company Limited, Japan and Professional Body Composition Analyzer IOI 353 of Jawon Medical, Korea.

### Body Composition Analyser Model IOI 353

**Make: Jawon Medical, Korea**

- Portable Body Composition Analyzer.
- 1 Machine = 32 Results
- Sensational Design, Reliability & Convenience.
- The Product uses 5 Factors, viz. Weight, Height, Gender, Age and Impedance to analyze the whole body.

### Uses:

- Best products for Brand Centers of MLM companies, Hospitals, Clinics, Gyms, Slimming Centers, Spas, Hotels, Resorts, BIO Medical Departments of hospitals etc.
- Best product for personal use, Doctors, Professionals, Trainers, Sports & Athletes, Patients, etc.

### Certificates & Website

Body Composition Analyzer IOI 353 is certified by CE, FDA, EN ISO 13485 & the reports are accredited with WHO and NIH standards.

For more details about us, you can visit us at [www.aarnasystems.in](http://www.aarnasystems.in), [www.indiamart.com/aarnasystems](http://www.indiamart.com/aarnasystems), [www.bonedensitometer.in](http://www.bonedensitometer.in), [www.aartech.co.in](http://www.aartech.co.in), and [www.aarnasystems.tradeindia.com](http://www.aarnasystems.tradeindia.com)

### Importer & Distributor :

## **AARNA SYSTEMS**

Hiran Magri, Sec.- 4, Udaipur-313 002 (Raj.) India

Phone : +91-294-2464136 Fax : +91-294-2464136

Cell : +9194603 28176, 097837 41286

Website : [www.aarnasystems.in](http://www.aarnasystems.in) | [www.indiamart.com/aarnasystems](http://www.indiamart.com/aarnasystems)

[www.bonedensitometer.in](http://www.bonedensitometer.in) | [www.aarnasystems.tradeindia.com](http://www.aarnasystems.tradeindia.com)

E-mail : [dkanthaliya@gmail.com](mailto:dkanthaliya@gmail.com) | [info@aarnasystems.in](mailto:info@aarnasystems.in)