

[ISO 13485: 2003, WHO-GMP & CE Compliance]

Office / Factory: A- 220, Street No: 5; Kabir Nagar; Delhi -10094(INDIA)

Phone No: + 91-11-22817003


Mobile No: + 91 – 11- 9555430720 & 9910178167

E-mail: info@medicalequipmentindia.com

Web: www.medicalequipmentindia.com

BROCHURE FOR NEURO SURGICAL PATTIES

**NEUROSURGICAL
PATTIES**



INTRODUCTION

Patties with thread are used to mop blood, fluids and as protective cover for delicate tissues in Neurosurgery and many other specialties. Purpose of thread is to identify the inserted patty.

MEI has developed world class patties. These are made of spun laced / thermally bonded non woven fabrics. Most of them have viscose as main fiber. These are highly absorbent.

Currently three types of base material and two types of thread are used. Thus there are six types of patties.

Material	Prefix	Example
Thin	No prefix	WM—
Thick	T Prefix	TWM—
Soft thin	P Prefix	PWM—

ADVANTAGES	Thread Example Suffix			
<p>1. Mechanically/thermally bonded fibers reduce shedding of fibers, reduced linting.</p> <p>2. Bonded fibers make removal of patties easy.</p> <p>3. Radio-opaque types can be identified in xray pictures.</p>	Black Nylon	No -WM —		
	Blue PP R-opaque	R -WM—R		
	PRESENTATION			
	<p>Neurosurgical patties are supplied in sterile peel open packs on a holding card for easy handling. Following sizes are regularly available, other sizes can be manufactured on demand.</p>			
DESCRIPTION				
<p>Non radio-opaque types have black nylon thread bonded with the patties. Radio-opaque types have 'blue polypropylene with bariun sulphate' thread bonded with the patties. MEI has developed a special process to bond the thread with patties. The Special automatic machine, developed in house by MEI, can bond and cut many lines of patties automatically.</p>	Code no.	Size (WxL)		
	Number in a pack	WM0505	5 x 5 mm	10 patties
	WM1010	10x10 mm	10 patties	
	WM1015	10x15 mm	10 patties	
	WM2020	20x20 mm	10 patties	
	WM2030	20x30 mm	10 patties	
	WM2040	20x40 mm	10 patties	
	WM2050	20x50 mm	10 patties	
WM5070	50x70 mm	3 patties		
WM70100	70x100 mm	2 patties		
	WARNING:			
	The Thread is for Identification only.			